- GB**
- USA**
- D**
- F**
- E**
- I**
- P**
- NL**

Instructions for use
Yasargil aneurysm clips

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




Legend

- A Test point
- B Force
- C Length of jaw part

Symbols on product

	Sterilization using radiation
	Not for reuse in intended applications as defined by the manufacturer
	Use by
	Attention, see instructions for use

Intended use

"Aesculap Yasargil aneurysm clips" will be referred to as "aneurysm clips" in the following text.

Special indication for permanent aneurysm clips only

The permanent aneurysm clips are intended for permanent stasis of cerebral aneurysms.

The permanent aneurysm clips are intended for single use only.

Contraindications: The permanent aneurysm clips are contraindicated for all applications except for the permanent stasis of cerebral aneurysms.

Special indication for temporary aneurysm clips only

The temporary aneurysm clips are intended for temporary stasis of cerebral aneurysms. Temporary aneurysm clips must not be implanted permanently under any circumstances.

The temporary aneurysm clips are intended for multiple use and can be reused repeatedly.

Contraindications: The temporary aneurysm clips are contraindicated for all applications except for the temporary stasis of cerebral aneurysms or blood vessels. Temporary aneurysm clips must not be used for permanent stasis or implantation.

Special Reuse Precautions for the US Market

Each aneurysm clip is individually packaged and provided sterile. In accordance with the Food and Drug Administration's (FDA) clearance requirements, each aneurysm clip is intended for use on a single patient during a single procedure. Temporary aneurysm clips are indicated for temporary occlusion of vessels and, therefore, repeated or multiple application during the same procedure will not adversely affect the efficacy of the aneurysm clip. Permanent aneurysm clips are indicated for permanent occlusion of aneurysms and, therefore, repeated or multiple application during the procedure should be avoided. Any aneurysm clip that has come in contact with a patient's blood or bodily fluid should not, under any circumstances, be cleaned, re-sterilized or in any other way prepared for use in another patient.

Product description

The aneurysm clips are available in two different materials:

- Cobalt alloy (Phynox) ISO 5832-7 and

• Titanium alloy Ti6Al4V, ISO 5832-3.

Titanium aneurysm clips and Phynox aneurysm clips are non-ferromagnetic. They are MR-proof to ISO 9713 in fields of up to 3.0 tesla, meaning the magnetic fields occurring during an MR tomography do not present any hazard. However, artifacts may occur. Such artifacts are more significant with Phynox clips than with titanium clips.

Stronger magnetic fields or an extension of the imaging field can lead to a dramatic increase in the incidence of positioning errors and artifacts. Depending on the MR pulse frequency, the size of the artifacts can vary considerably.

Different models of our aneurysm clips are available. Please contact Aesculap for further information and ordering of special brochures about titanium and Phynox aneurysm clips.

The aneurysm clips are available in three sizes (Mini, Standard and Long).

The closing force of each aneurysm clip is measured individually and printed on the individual packaging. The closing force given on the packaging was measured at 1/3 along the length (from the tip) of the jaw part, at the center of the contact area. The Phynox aneurysm clips are measured with the jaws opened to 0.5 mm (see Fig. 1) whereas the titanium aneurysm clips are measured with the jaws opened to 1 mm (see Fig. 2).

Each aneurysm clip carries an individual serial number.

The aneurysm clips are color-coded to differentiate the various sizes and applications.

Special product description for permanent aneurysm clips only

	Phynox aneurysm clips			Titanium aneurysm clips	
Size	Mini	Standard	Long	Mini	Standard
Color coding	silver			Jaw part: silver Spring: pink	Jaw part: silver Spring: blue

Each permanent aneurysm clips package contains additional labels showing the article number and the individual serial number of the clip.

Special product description for temporary aneurysm clips

	Phynox aneurysm clips			Titanium aneurysm clips	
Size	Mini	Standard		Mini	Standard
Color coding	gold			Jaw part: gold Spring: pink	Jaw part: gold Spring: blue

- Please refer to special reuse precautions for the US market.

Safe handling and preparation

CAUTION

Federal law restricts this device to sale by or on order of a physician!

- Read and observe the instructions for use, and keep them in a safe place.
- The product may only be used for its intended application, see Intended use.
- Prior to each use, inspect the product for: loose, bent, broken, cracked, worn, or fractured components.
- Do not use the product if it is damaged or defective. Products found to be damaged must be set aside immediately.

Each aneurysm clip is individually packaged and sterilized by radiation (min. dose 25 kGy).

If the sterile wrapping is open, torn, perforated or otherwise damaged, or when the sterility of the aneurysm clip has "expired" (use-by date)", the aneurysm clip must be regarded as unsterile and has to be sterilized according to the instructions given in the sections on cleaning and sterilizing.

Each individually packaged aneurysm clip is supplied in double sterile packaging including instructions for use and labels.

Special handling and preparation for permanent aneurysm clips

Aesculap recommends preparing the permanent aneurysm clips prior to operation in their unopened sterile packaging. In this way, damage to the booster clips is avoided, their functionality is ensured and their correct closing force is maintained.

Labels

Each permanent aneurysm clips package contains additional labels showing the article number and the individual serial number of the clip. To facilitate the subsequent radiological examination, the patient's file (for the hospital) and the patient ID (for the patient) should be marked with these labels.

Patient ID

Important information concerning the aneurysm clip implanted and the operation can be noted down on the patient ID. To facilitate the postoperative radiological examination, each patient should receive a patient ID.

- Contents of the patient ID
- Patient data

– Hospital data

– Date of operation

– Name of surgeon

– Article number and individual serial number of the aneurysm clip implanted

The patient ID can be ordered separately from Aesculap at a later date.

To exclude damage, inaccurate functioning and incorrect closing force:

➤ Apply each aneurysm clip with the appropriate applier and removal forceps.

➤ Never use the aneurysm clips with applier or removal forceps from other manufacturers

- The Phynox aneurysm clips must be applied with the applier for Phynox aneurysm clips. The titanium aneurysm clips must be applied with the applier and removal forceps for titanium aneurysm clips. Mini, Standard or Long aneurysm clips may be used only with the applier and removal forceps of the appropriate size (Mini, Standard or Long).
- All applier and removal forceps are marked according to size (Mini, Standard or Long) and clip material so that their correct application with clips of the respective size and material is ensured. Additionally, the applier and removal forceps for titanium aneurysm clips can also be color-coded according to clip size.

For further information on the suitable applier and removal forceps, please contact Aesculap, which can also supply you with special brochures for titanium and Phynox aneurysm clips.

- Carefully pick up the whole aneurysm clip with the jaw piece of the appropriate applier forceps (see Fig. 3). The aneurysm clip can be damaged, jump out or slip if picked up improperly (see Fig. 4/5/6).

To avoid possible risks to the patient:

➤ The permanent aneurysm clips are designed for single use only.

➤ Temporary aneurysm clips must not be implanted permanently under any circumstances.

Aesculap cannot accept any responsibility for aneurysm clips that are handled inappropriately or not according to the present instructions for use.

Safe handling

Aneurysm clips should only be applied by surgeons appropriately trained for, and experienced in, the use of aneurysm clips.

If an aneurysm clip appears changed or shows signs of damage, e.g. incorrect jaw position, bent parts or changed closing force:

➤ Set aside the aneurysm clip.

To avoid damage to the aneurysm clips:

- Always handle the aneurysm clips with appropriate care.
- Never open an aneurysm clip with your fingers.
- Avoid manual and/or mechanical manipulation of the aneurysm clip.

Excessive, rough or repeated handling, especially opening and closing of aneurysm clips, be it in general use or during cleaning and sterilization, can change the closing force and impair the clinical effectiveness of the aneurysm clips.

- As soon as the sterile aneurysm clip has been brought to the operation field in the usual manner, pick up the aneurysm clip from its double sterile packaging and position it carefully between the jaws of the applier forceps.

The correct positioning of the aneurysm clip between the jaws of the appropriate applier forceps (see Fig. 3) is crucial here. If the aneurysm clip is not placed precisely between the jaws parts of the appropriate applier forceps (see Fig. 4/5/6), this can lead to damage to the aneurysm clip. It could cause a change in the factory-set closing force or result in the aneurysm clip slipping from the applicator, which could be dangerous during an operation. These directions for use must also be observed when using the removal forceps.

- To preclude the possibility of a galvanic reaction, avoid implanting permanent aneurysm clips made of different materials and/or by different manufacturers in such a way that they touch each other.
- Make absolutely certain that the aneurysm clip sits correctly on the aneurysm neck and close to the blood vessel, during and after implantation.

The application of aneurysm clips involves the following severe risks:

➤ Shifting or breakage of the aneurysm clip

➤ Yawing of the jaw parts

➤ Rupture of the aneurysm due to punctual, incomplete contact of the clip jaws on the aneurysm neck.

➤ For large aneurysms, reduction of the blood vessel cross-section close to the edge of the blood vessel.

➤ Cerebrovascular spasms and sudden death.

➤ Infections of the operation wound and general surgical complications are other unwelcome (side)effects.

Each patient must be informed comprehensively about the properties of aneurysm clips, as well as about the surgical procedure.

Care and handling

Each aneurysm clip in its unopened original packaging has been packed and sterilized individually and is supplied as a sterile product. Unsterile aneurysm clips must be sterilized following the cleaning and sterilization processes outlined below.

To preclude damage, inaccurate functioning and an incorrect closing force:

➤ Use each aneurysm clip with the appropriate applier and removal forceps.

➤ Never use the aneurysm clips with applier or removal forceps from other manufacturers.

Excessive, rough or repeated handling, especially opening and closing of aneurysm clips, be it in general use or during cleaning and sterilization, can damage the aneurysm clip, change the closing force and impair the clinical effectiveness of the aneurysm clips.

To avoid damage to the aneurysm clips:

- Always handle the aneurysm clips with appropriate care.
- Never open an aneurysm clip with your fingers.
- Avoid manual and/or mechanical manipulation of the aneurysm clip.

Cleaning and disinfection

Special cleaning instructions for permanent aneurysm clips

Permanent aneurysm clips that have been implanted before or been in contact with blood, tissue, or body fluids:

➤ Must be set aside for disposal.

➤ Must not be sterilized.

Special cleaning instructions for temporary aneurysm clips

Note

For patients with Creutzfeldt-Jakob Disease (CJD), suspected CJD or possible variants of CJD, the relevant national regulations concerning the reprocessing of the products must be observed.

If necessary, sterilization must be preceded by the following procedure:

- Carry out thorough and careful manual pre-cleaning of the aneurysm clip, with or without ultrasound treatment.
- Use an appropriate applier forceps to open the jaws parts of the aneurysm clip.
- Carry out a final rinse with distilled or demineralized water, making certain the all remaining debris was removed.

Inspection

- Inspect each individual aneurysm clip.
- Any aneurysm clip found with one of the following characteristics must be set aside and excluded from further use:
 - Signs of damage
 - Incorrect jaw position
 - Bent components
 - Changed closing force
 - Misalignment
 - Soilage that cannot be removed

The color coding may fade in the course of processing. Aneurysm clips whose color coding cannot be identified unambiguously anymore must be set aside and excluded from further use.

Sterilization method and parameters

Each aneurysm clip in its unopened original packaging has been packed and sterilized individually and is supplied as a sterile product.

CAUTION

An aneurysm clip booster clip must be regarded as unsterile e.g. if the sterile packaging is open, torn, perforated or damaged, or if the sterility of the product has "expired" (use-by date).

- Sterilization has to be carried out through a validated sterilization process.
- Using a non-validated sterilization process may result in the aneurysm clip being unsterile.
- Sterilize the aneurysm clip in the appropriate storage devices recommended by Aesculap for Aesculap Yasargil aneurysm clips. For further information on the appropriate storage devices recommended by Aesculap, please contact Aesculap or consult the special brochures for titanium and Phynox aneurysm clips. These brochures too can be ordered from Aesculap.

➤ The sterilization must be carried out according to a validated steam sterilization process, e.g. in a sterilizer conforming to standards EN 285/ANSI/AAMI/ISO 11134-1993, ANSI/AAMI ST46-1993 and validated according to EN 554/ISO 13683. For the fractionated vacuum process, the following parameters must be observed: 134 °C (2-bar program) with a minimum holding time of 5 minutes.

Sterilization for the US market

- Sterilization of the device may be accomplished by steam.
- Aesculap does not recommend the device be sterilized by "Flash" or chemical sterilization.
- Surgical instruments may also be placed within an Aesculap rigid sterilization container (sterile container) for processing under generally accepted hospital in-use conditions.

The recommended sterilization parameters are as follows:

Sterilization method	Temp.	Minimum exposure time	
		Wrapped	In a sterile container system
Pre-vacuum	270—275 °F	5 min	5 min

Storage

Store the aneurysm clips:

- either in their sterile packaging or
- in a compatible storage system recommended by Aesculap, kept in a dust-free area protected from chemical fumes and extreme air or temperature fluctuations.

For further information on the appropriate storage devices recommended by Aesculap, please contact Aesculap or consult the special brochures for titanium and Phynox aneurysm clips. These brochures too can be ordered from Aesculap.

Distributor in the US/Contact in Canada for product information and complaints

Aesculap Inc.
3773 Corporate Parkway
Center Valley, PA 18034
USA

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